


How to make the low-tech WordPower 108 w/Keyboard book with flips



Why use a low-tech option?

- * Use as a **back - up** system for a high- tech system.
- * Use as a **beginner** board. This allows you to transition easily to the high- tech option using the same vocabulary arrangement. Or used while waiting for high-tech device.
- * Use as a **model** for teaching a classroom of students with language deficits.
- * The **communication partner** can also model using the low - tech option.
- * Use during **water** activities. If you punch holes in it, probably not great for submerging, but would last a decent amount of time to water splashes (water table). If you wanted to make the hole punched sheets more waterproof, leave enough lamination above the sheets to punch holes.


What materials are recommended?

- ✓ **PDF file** of the book
- ✓ **Printer**-preferably color
- ✓ **Lamination** – For the book pictured here, I used a home laminator with 3 ml heat sealed sheets. I have used this laminator for years (over 10) for many different things (placemats, swimming games charts, etc.)
- ✓ **Binding** – Use 3 individual rings or a binder. I chose to leave the binder as is and added a felt stand to the outside of the book (attached with Velcro). This allowed the flips to stand up and then were easily flipped down instead of having to flip from the back. You could select any binding option that suits you. Basically you will want the flips to be able to flip up.
- ✓ **Hole punch**- a 3 hole punch works nicely. This will be determined by the method of binding.
- ✓ **Scissors**



How?

1. **Print** out PDF file.
2. **Add** any **custom** pictures or words to appropriate flip locations. You can quickly glue a small square on a blank or over a word you do not find necessary for your client. It will be laminated over and kept secure.
3. **Cut** out the pieces. The full page vocabulary should sit at the correct position to align the flips the way they are placed on the pages. You will see a dotted line between the flips to guide you in cutting. This cutting guide will provide you with the top of the next flip. Cut around the flip tabs and follow the line over, as well as straight across the bottom of the other flips. This should align the tops of the flips and leave equal space for the binding option. Printers vary in the way they align and print so check to see that it aligns correctly when you begin cutting. Note: the ACTIONS page will be a whole page flip that sits behind the main page. Cut the tab on main page or straight along bottom.
4. **Laminate**- If you are using a page lamination system, you will simply put the pieces in the pages allowing enough space to leave a laminated edge. If you are using a roller type of laminator, make sure you leave enough space around the edges for the cut and a laminated edge.
5. **Cut** the pieces again- leave a small amount a lamination around edges to keep a seal (lasts longer this way)
6. **Bind**-Use a 3-hole punch or other punch option. (Rings or binder) Check to see the hole alignment matches as you go. Or use a binding machine.
7. **Other** -To make the binder stand up- you can add a strap or piece of material to back and front of the folder with Velcro.



Word power

Based on WordPower 108 w/keyboard



PEOPLE 	QUESTN 	PLACES 	TIME 	SOCIAL 	DESCRB 	good 	more 	-ed	-ing	-s	GROUPS
I	me	Q	W	E	R	T	Y	U	I	O	P
my	don't	A	S	D	F	G	H	J	K	L	bksp
it	not	Shift 123	Z	X	C	space	V	B	N	M	.
he	is	to	be-	come 	eat 	drink 	feel 	at	and	a-	the
she	are	do	did	go 	get 	help 	like 	for	in	on	that
they	was	can	could	make 	need 	play 	put 	read 	of	off	this
we	have	will	would	say 	see 	stop 	take 	talk 	up 	out	any- or every- if
you	your	ACTIONS 	want 	tell 	think 	turn 	watch 	work 	down 	with	clear

Cut out area for tab or remove entire bottom white area




people 	him	his	family 	mom 	dad 	grandma 	grandpa 	aunt 	uncle 	parent 	person 
mine	her	hers	friend 	them	their	sister 	brother 	girl 	boy 	woman 	man 

1


PEOPLE


										teacher 	student 	

2

has	?	had	were		what	when	where	SOCIAL QUESTNS 		question 
should	may	does	why	how	which	who	whose	Guess what?		What will we do... 

3

QUESTN


building 	airport 	aquarium 	bank 	bowling 	church 	doctor 	gas stat 	grocery 	hospital 	hotel 	place 
		library 	mall 	movie 	post off 	store 	rest'rant 				

4

PLACES


home 	house 	bathrm 	bed 	bedrm 	closet 	dining 	garage 	kitchen 	laundry 		
		living rm 	office 	room 	window 						

5

outside 	beach 	camp 	farm 	garden 	island 	lake 	ocean 	playgrnd 	park 		GEOGRA 
		pool 	woods 	yard 	zoo 						


6

school 	bus 	class 	library 	lunch 	therapy 	speech 	OT	PT			
					baseball 	basketball 	football 	soccer 	tennis 		

7

Sun	Mon	Tues	Wed	Thurs	Fri	Sat	last 	next 	yest'day 	today 	tom'row 
time 	morning 	aft'noon 	day 	night 	week-end	week	before 	after 	now 	soon 	later 
			TIME 								

8

Jan	Feb	March	April	May	June	July	Aug	Sept	Oct	Nov	Dec
month		2016	2017	year			season 	spring 	summer 	fall 	winter 

9



10



























11



12



broken 	fixed 	light 	dark 	long 	short 	hungry 	thirsty 	tired 	sick 	fun 	right 
clean 	dirty 	cold 	hot 	loud 	quiet 	excited 	nice 	proud 	same 	diff'nt 	wrong 

13

cool 	warm 	near 	far 	new 	angry 	bored 	frustr'd 	mad 	mean 	pretty 	ugly 
fast 	slow 	full 	empty 	old 	scared 	uncomf 	wet 	dry 	smart 	dumb 	

14

best 	worst 	few 	many 	early 	late 	excellent 	terrible 	fat 	thin 	most 	least 
cheap 	expensive 	first 	last 	high 	low 	exciting 	boring 	more 	less 	true 	false 

15

afraid 	alone 	depress 	embarrass 	lonely 	nervous 	surprised 	terrific 	worried 			
apart 	together 	young 	old 								

16

color 	red 	orange 	yellow 	green 	blue 	purple 	silver 	round 	square 	triangular 	rectangul 
	black 	white 	pink 	brown 	grey 	tan 	gold 	diamond 	octagonal 	straight 	crooked 

























17

DRINKS 	water 	juice 	milk 	choc mlk 	soda 	iced tea 	lemonade 	shake 	smoothie 		
SNACKS 	applesau 	crackers 	candy 	cookie 	fruit snack 	peanut b 	popcorn 	pretzels 	goldfish 	yogurt 	chips 
											GROUPS 

18

MEALS 	sandwich 	macaroni 	cheese 	fr fries 	hamburg 	hot dog 	nuggets 	pizza 	spaghet 	salad 	soup 
cereal 	oatmeal 	toast 	eggs 	pancake 	syrup 	peanut b 	jelly 	bacon 	sausage 	grilled ch 	chicken 

























19

fruit 	apple 	banana 	grapes 	blueberr 	orange 	cherries 	grapefru 	raspber 	strawber 	pear 	waterme 
veggie 	broccoli 	corn 	gr beans 	peas 	carrot 	celery 	cucumb 	lettuce 	tomato 	potato 	pepper 

20

art 	crayon 	marker 	pen 	pencil 	paintbrsh 	paint 	glue 	scissors 	tape 	paper 	picture 
pets 	bird 	cat 	dog 	fish 							


21

body 	head 	ear 	eye 	nose 	mouth 	arm 	leg 	hand 	chest 	back 	stomach 
health 	sick 	cold 	fever 	headache 	sore throat 	stomach 	toothache 	medicine 	doctor 	dentist 	nurse 

22



23



24



25

if	then	again	all	around	away	front	back	behind	top	above	over
here	there	these	those	up	down	left	right		bottom	below	under

GROUPS

P

bksp

.

the

that

this

any- or
every- if

