
TWIN POWERS UNITE READING PROJECT! Teaching AAC Core AND Dolch Sight Words



Housekeeping Info



ASHA Disclosures

- Maureen Castillo, M.A., CCC-SLP
 - Financial-none
 - Non-Financial ASHA Member
- Lisa Timm, M.S., CCC-SLP
 - Financial- employee-owner of Saltillo Corporation
 - Non-Financial ASHA Member, ASHA SIG 12, Iowa Speech and Hearing, USSAAC Board and member



Sable Gegeneties is approved by the Contracting Education Record of the American SpaceUninguege/Hearing Association (RGM) to provide certificity advoction and other another of ASMA (Dis), indexeducing thes associal information for method of ASMA (Dis), indexeducing their associal information for method of ASMA (Dis), indexeducing the advocation account, specific products or distant procedures.

This course is offered for .10 ASHA CEUS. (Beginner level, Professional area)

ASHA Housekeeping Info...

- Email with post-test within one hour
- Submit within 15 days
- Must stay online the entire time
- Instructions and participation forms are in the materials section

Who Am I???

Maureen Castillo, M.A., CCC-SLP

maureen.castillo@adams12.org

-SLP for 27 years, 17 specializing on AAC in the public schools

-Assistive Technology Specialist Adams 12 Five Star Schools

-Owner/Provider of iTherapy, LLC

Implementing an AAC Device in the Schools

What are the Challenges?????

Challenges...


...are what makes life interesting; overcoming them is what makes life meaningful.

- Josina J. Herber

Challenge 1: Where do I start? What is the focus?



Challenge 2: Training Multiple People



Challenge 3: Generalizing ! From Therapy to Real Life!


In School ! At home! With peers!




Challenge 4: Including AAC in Literacy

Access to books our students can independently

read or share in the reading is a challenge!





CORE Language (Pillar one)

- Small set of commonly used words that support communication and language learning
- 85% of what we say is communicated with only 200 basic core words
- Mostly pronouns, verbs, descriptors and prepositions very few nouns!
- Consistent across place, topic, and cognitive ability

- For a child to become an effective communicator, either verbally or with the assistance of a device, he/she must master these basic core words

DLM First Forty

- Dynamic Learning Maps (DLM) and the Alternative Assessment Consortium based on work by Karen Erickson at the University of North Carolina.

- These words are recommended for students who use augmentative and alternative communication (AAC)

 These words were selected from a larger core vocabulary set based on their usefulness in everyday communication and in addressing the DLM[™] Essential Elements.

DLM'sTM "First 40"

1.	I	like	not	want	
2.	help	it	more	different	
3.	who	she	you	he	
4.	where	up	on	in	
5. me		make	get	look	
6.	what	need	are	is	
7.	some	put	all	this	
8.	don't	that	go	do	
9.	when	finished	can	here	
10.	open	turn	stop	over	



DLM's First 40

- Previous slide shows our first 40 suggested words divided into 10 groups of 4 words each.
- The vocabulary selection IS research-based, however, the organization is NOT.
- The grouping of 4 words is one approach to organization.
- Organization was based on providing multiple opportunities for expression as well as receptive input and modeling IF a student could handle only 4 symbols at a time.

Center for Literacy & Disability Studies

University of North Carolina at Chapel Hill



Vocabulary Decisions: **Continuum of Usefulness**

Useful for restricted Useful for a wide purposes in limited range of purposes in contexts on a a variety of contexts specific topic on numerous topics isosceles trapezoid measure socialize dog Pittsburgh Zoo forest Harley want star nibble Core up opinion Fringe restaurant you angle meteorite banana peppers opposite shoes summarize more turn protect park wonder march **Red Hot Chili Peppers** circle proud explore

Lots of opportunities to teach throughout the day.

Opportunities limited to specific activities or specials.



Center for Literacy & Disability Studies

University of North Carolina at Chanel Hill

Dolch Sight Words (Pillar Two)


- Dolch Word List, was compiled by Edward William Dolch, PhD, in 1948
- The list was originally published in his book "Problems in Reading"
- 220 "service words" that must be quickly recognized in order to achieve reading fluency.
- The Dolch Word List is also called Sight Words or The Dolch 220.
- It includes the most frequently used words in the English language.
- Sight words make up 50 to 70 percent of any general text

Which is more Important???

Matter of perspective.....

Communication - CORE

Reading - Dolch/Sight Words



Valuable Real Estate

If you could only speak 40 words, would you want one of them to be "the"?



Twin Powers Unite - Working Together!

Twenty Words We Have in Common - Cross Referenced

(20 on DLM First 40 list AND 20 on Dolch Pre-Primer or Primer Lists)

Chosen Based on Scholastic Sight Word Readers

Although 20 Words are Targets, due to other high frequency words in texts, many more will be covered (ie. see, play, get)

This is where you start!!!



1	The second	like	not	want
2	help	it	more	different
3	who	she	you	he
4	where	up	on	in
5	me	make	get	look
6	what	need	are	is
7	some	put	all	this
8	don't	that	go	do
9	when	finished	can	here

and the second	Dai	mer	Grade One	Grade	Two	Grade Three
Pre-primer	Pro	mer	Crass	A PERSONAL S		about
ALL CONTRACTOR	all	under	after	always	why	better
a and	am	want	again	around	wish	bring
	are	was	an	because	work	
away big	at	well	any	been	would	clean
blue	ate	went	ask	before	write	ciedri
can	be	what	as	best	your	done
come	black	white	by	both		draw
down	brown	who	could	buy		drink
find	but	will	every	call		
for	came	with	fly	cold		eight
funny	did	yes	from	does		fall
go	do	1	give	don't		far
help	eat	1942-64	going	fast		full
here	four	1. 1. 17 1. 19	had	first		got
I	get		has	five		grow
in	good		her	found		hold
is	have		him	gave		hot
it	he	had a	his	goes		hurt
jump	into		how	green		if
little	like		just	its		keep
look	must		know	made		kind
make	new		let	many		laugh
me	no		live	off		light
	now		may	or		long
my not	on		of	pull		much
one	our		old	read		myself
	out		once	right		never
play red	please		open	sing		only
1000	pretty		over	sit		own
run			05725752	0.20078		1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1
said	ran		put	sleep		pick
see	ride		round	tell		seven
the	saw		some	their		shall*
three	say		stop	these		show
to	she		take	those		six
two	50		thank	upon		small
up	soon		them	us		start
we	that		then	use		ten
where	there		think	very		today
yellow	they		walk	wash		togethe
you	this		were	which		try
	too		when			warm

Here were our Twenty Words! (on both lists)

1.	Ι	11. in
2.	Like	12. here
3.	Want	13. are
4.	lt	14. can
5.	Не	15. go
6.	You	16. help
7.	Up	17. get
8.	On	18. that
9.	Look	19. is
10	. Me	20. this

Partner Augmented Input (Modeling) - Pillar Three

- A powerful, evidence based approach in which communication partners point to the symbols on the child's communication device while simultaneously talking
- As each book is introduced, an adult will "read" book with the AAC device following a guide created by Chat Editor
- Peers will also use this guide to learn and read using their friend's talker!
- How we teach an "umbrella full of people" the device!

Video on Using Modeling

Power of Peer Modeling - Pillar Four

- Strategy researched by Kathy Bourque from the University of Kansas
- Children responded better to peer modeling than they did to adults!
- Reading partners or book buddies as peers
- High expectations for peer to be in role as "teacher"
- Goal is to improve sight word recognition of both students
- Using book as common denominator, helps us generalize




Teaching and Modeling Vocabulary

- No student will be able to use this approach without teaching and modeling.
- Teaching:
 - Teach the words in the set you select.
 - Teach the symbols in the set you select.
- Modeling
 - Model the use of the vocabulary during instructional activities and all other communication activities.
 - Receptive input using the system is critical.
 - Create multiple systems and take advantage of peers.

Center for Literacy & Disability Studies

University of North



Anchor Books - Pillar Five

Chose a commercial set of books to create


paths for as a basis for the target words each

week. Liked that students had a "real"

book, just like their peers!

Why? Fewer than 10% of AAC users read

beyond the second grade level (Erickson, 2003)



Who was in our Project?

- Ten students targeted in grades K-1 who use either NovaChat or TouchChat to communicate (20,42 or 60 layout) or PRC device with 45 Sequence
- Each student was paired with a peer who also needs to learn his/her sight words (RTI students)



Data Collection

-Baseline Data will be collected on both peers

- Three data points to collect:
 - 1. Receptive knowledge of targeted words in field of 3
 - 2. Receptive knowledge of targeted words in phrase or sentence
 - 3. Expressive (verbal or on device) recitation of targeted words

Data sheets provided in packet. (Receptive, Context and Expressive)



Reading from Field of Three


Reading from a Field of

Three Video Link



Reading in Contextual Sentence

Sentence Video



Expressive - Reading with Voice or Device





Project Schedule

Two books/week with 2-3 targeted words - 5 min/day - schedule in packet


Each child receives copy of each book

Adult reads selection to both students using provided guide on Monday

Two times/week adult reads again with AAC user - may be shared reading where student fills in target word

Two times/week both peers read selections together

Weekly intervention data sheet



Project Implementation Steps

- 1. Choose target words and books
- 2. Choose students and peers
- 3. Collect baseline data
- 4. Create paths for books
- 5. Create timeline
- 6. Implement!
- 7. Collect post intervention data
- 8. Analysis



Schedule Targeted Words for AAC and Sight Word Project:

DLM "First 40" Core Word AND Dolch Pre-Primer/Primer Words (on both lists):

Word Take Pre-Data/	<u>Book</u> Training	<u>Week</u> January 11-15	Theme/Concept
	Farm Friends	1 - January 19-22	Animals/Actions
like	Animals		
want	New Socks	2 - January 25-29	Clothing/Actions
it	He Did It		
he			
you	My Dog	3 - Feb. 1-5	Actions/Positions
up	Up and Down		(Rhyming)
		Review	
on	Ride On!	4 - Feb. 16-19	Vehicles/Colors
look	City Colors		
me	Me Too!	5 - Feb. 22-26	Actions/Insects
in	Guessing Game		(Inferencing)
here	The Party	6 - Feb. 29-March 4	Party
are	• •		
can	My Bear	7 - March 7-11	Actions/Places
go	Go, Go		
help	Helpers	8 - March 14-18	Places/Food
get	Dinner		
that	That Hat	9 - March 21-25	Describing/Groups
is		and the second second	(Possession)
this	Little and Big This Is A Peach		(*
Take Post Data		March 28-April 1	
20 Words	18 Books	9 instructional weeks of	is 2 data collection weeks

Reading Paths and Guides

Reading paths for each book will be provided on a separate sheet to promote independence

Reading guide supplied to block out


extraneous information


At end of project, paths with guides

can be sent home for parents to use

Might be used in ESY as well!

Helps us model and generalize!






A note about the paths







Visual Supports using Chat Editor

This is a FREE download from Saltillo website : www.saltillo.com


What is **Chat Editor**?

Chat Editor is a supplemental program that has been provided for customizing the NOVA chat vocabulary on your Windows desktop or laptop computer. Although customizing is also possible on the NOVA chat device itself, the Editor allows the family and/or professionals to customize the vocabulary when the NOVA chat device is not present.

Note: Having the custom vocabulary on a computer is always recommended as a backup to the NOVA chat device

Chat Editor

Select the vocab you want to open.





Your vocab is then displayed.

Chat Editor (42 Basic)									
Settings L	Settings Library Help								
	🚱 🚱 😭 🧶 Edit Mode Capture								
ABC [►] 123									
clear	I	me	don't	to	be-	come Come		a	more
•	my	can 🗅	eat D	9		go 🗅	this	the D	that
delete wd	it	do	get 🎴	h sr	know P	•	about	and	at
-5	you	is D	like D			STOP	for	in D	with
we	your	∎ have	think P	want 🍳	watch ^D	work P	of	on	out


Select Capture to create you icon sequence. Click on the buttons to create you sequence.

Chat Editor (42 Basic)									
Settings Library Help									
	C C T C Capture								
Settings I	Library He	lp							
3	1	X Edit	Mode Cap	ture					
1									
ABC [►] 123					buy D	call □	really	got 🗅	went ⁱ
	would	will	don't	to	be-	come D	(⁽²⁾ ©		÷,
should	could	can	eat P	<u></u>		go 🏼	*0	look D	
delete wd	did	do	get D	help ▷	know 🍳	listen 🍳	make P	put &	
-'m	was	am	like D	need D	play D	stop P		ि	see I
ACTIONS A - Z	had	have	think D	want P	watch ^D	work	take D	talk D	tell ^I
Button Capture									
Ι									
	Fyit Copy To Clipboard Clear Clarture Selections								





Repeat the process of copy to Clipboard and paste to Google Slide



When finished; you can Print or Share your strips..

Reading with my Teacher

New Socks Book



With my Paraeducator

Paraeducator reading with student



With a Peer!

Partner Modeling



In the Classroom!

In the Classroom




Commenting

At the end of each week, both students

were asked to complete a comment

sheet for each book.



Why? Commenting provides each student with a unique voice and way to state their opinion! This is often a forgotten communicative function for our AAC users!

Post Intervention Data Collection

Data was collected at end of 9 weeks - really 8 due to snow days!

Students took all books home with reading guides



Book Bag

Each child was provided with a book bag to store books, flashcards and comment sheets.

At end of project, books were sent home and parents encouraged to read with their child! Again, modeling and generalization!



Data Collection Tools

-Project Information Summary Sheet

Pre/Post Data

Demographic Information

-Target Words/Data Collection Sheets

Expressive

Receptive


Context (Sentences are provided)

-Weekly Tracking Sheet

Used by all groups using the program

Place for comments/extra data collection

Our Results




7 out of 10 students showed growth, with 3 students with mastery of receptive sight words both pre/post.


Receptive Word Acquisition

Range = 0-13 words learned

Average = 4.7 words learned



Our Results




9 out of 10 students showed growth in expressive labeling of target words using either oral output or AAC device or signing, 1 student showed mastery both pre/post.

Expressive Word Acquisition


Range = 0-10 words learned

Average words learned = 4.8

2nd cycle Pre-expressive and Post-expressive




Our Results



9 out of 10 students showed growth in identifying words within context (simple sentences), 1 student showed mastery both pre/post.

Words in Context

Range = 0-13 words learned Average = 5 words learned




Overview of Peer Comparison

Peers showed either growth or maintenance of skills across areas. 5 out of 6 peers were higher or comparable to student in project. In one peer/student pair, student showed greater growth than peer.

Peers provided great models in some cases, however in others gained needed literacy skills.

POST Twin Powers Unite Survey Results

How Helpful was ...



Top Three Take Aways from Staff

- 1. Twin Powers Project was an aid for device generalization!
- 2. Modeling was an effective tool for staff to learn the device!
- 3. Providing the paths allowed staff to model!



Feedback - What worked?

ALL students improved!!

Modeling as a tool to include paraeducators, teachers, peers, worked!

Anchoring core to a book helped create literacy opportunities!

Project as a guide for generalization of device usage!

Overwhelmingly, everyone wants program to continue/expand!

Feedback - What would we do differently?

Pacing - one book/week - plan for snow, assessment days, illness, etc.

Group based on themes and suggest supplemental text to promote generalization language goals/themes

Sentence strips to work on sequencing sentences and writing extension ideas

Choose text which is user friendly! Writing these!

Core Word Wall

DLM-Word Wall Link



Bonus Ideas....

Classic Book - generalize targeted word recognition in print

Themes - build vocabulary based on categories and extend to writing

Units of Study - Grouped to help you incorporate lessons throughout day

Writing - cut apart sentences, worksheets, etc.

Bitsboard and other apps!


	Core Word of	Scholastic	Classic Book	Extension Book	Theme	Unit of Study
Dates	the Week	Reader				
Aug. 31	Training					
Sept. 6-9	Data Collection					
Sept. 12-16	me	Me Too		Who Can?	Actions	Life Science-Animal s
Sept. 19-23	like	Animals	Brown Bear by Eric Carle	What I Like!	Actions	
Sept. 26-30	1	Farm Friends	и и	How Do I Feel?	Animals	
Oct. 3-7	look	City Colors	и и	Look!	Colors	
Oct. 10-12	Review Week					
Oct. 17-21	go	Go, Go	Go Away Big Green Monster By Ed <u>Emberly</u>	Where Can You Go?	Places	Civics - Community
Oct. 24-28	up	Up and Down		What is Up?	Concepts up/down	
Oct.31-Nov. 4	want	New Socks		Surprise		
Nov. 7-10	can	My Bear	From Head to Toe - Eric Carle	We Like to Help	Body parts	Life Science – Human Body
Nov. 14-18	you	My Dog		Where?		
Nov.28-Dec. 2	it	He Did It	и и	It is Tie	Clothing	
Dec. 5-9	Data					

Extension Ideas - pair with text for Generalization

Brown Bear




Building Sentences



Predictable Chart Writing

9 2. ee a See 9 20 See See N



me	me	me		
O	O	O		
me	me	me		
O	O	O		
me	me	me		
O	O	⊕		

Our classroom word of the week


is



Please model the way and use our word of the week in your conversations with our students.

Thank You!

Cut Apart Sentences Farm Friends-WP42B



I	see	a	horse
---	-----	---	-------

Extension Books

Original books based

on core! Compliment


to instruction to help

differentiate as well as

great for secondary

students!


How Do I Feel?



Freebie for the Day

I am giving you all a free pdf copy of my original book, "How Do I Feel" as well as the accompanying reading guides for Word Power Simply 20, 42 Basic and 60 Basic. Enjoy!

How Do I Feel Original Book and Reading Guides



Thanks to All the Staff who helped with this Project!



Thank you for joining us!



Join us on Social Media!

