WordPower48 Español – teaching tips and support

Nancy Inman & Pamela Corley

1. Access to Quick Basic Communication with Picture Support:

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Example 2. me gusta "it pleases me" (similar to "I like" in English)

me gustaría "it would please me" (notice link to COSAS FAVORITAS "favorite things")

Me gusta ir con usted. "I like to go with you."	Me gusta (ACCIONES) escuchar música. "I like to listen to music."
Me gusta comer en Chick-fil-A y Taco Bell. "I like to eat at Chick-fil-A and Taco Bell."	Me gusta jugar con usted. "I like to play with you."
Me gusta mi (GENTE) amigo. "I like my friend."	Me gustaría (ACCIONES) leer un libro. "I would like to read a book."

Example 3. tengo "I have" (calor "heat" sed "thirst" frio "cold" hambre "hunger" miedo "fear" sueño "sleepy")

Tengo hambre "I am hungry." (literally "I have hunger.")	Tengo sed. "I am thirsty." (literally "I have thirst.")
Tengo sueño. "I am sleepy."	Tengo miedo. "I am afraid." (literally "I have fear.")

Example 4. me siento "I feel" (literally "I feel myself")

(link to DESCRIPCIÓN). Notice "m" for masculine and "f" for feminine

Me siento feliz. "I feel happy."	Me siento cansado (masc). "I feel tired."	Me siento cansada (fem).
-------------------------------------	--	--------------------------

Example 5. FRASES (quick access to often-used phrases).

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2. Direct Link from Verbs to Related Categories:

Example a) **comer** “*to eat*” (link to **COMIDA** “*food*”)

Quiero comer pizza con pepperoni y cebolla. “I want to eat pizza with pepperoni and onion.”	Quiero comer papas fritas con salsa. “I want to eat chips with salsa.”
Me gusta comer huevos con tocino. “I like to eat eggs with bacon.”	Quisiera comer con usted. “I would like to eat with you.”

Example b) **beber** “*to drink*” (link to **BEBIDAS** “*drinks*”)

“Quiero beber Coca-Cola y el café con crema.” “I want to drink Coca-Cola and coffee with cream.”	Quiero beber un poco de limonada. “I want to drink a little lemonade.”
Me gusta beber <u>refresco</u> de banana. “I like to drink banana soda.”	Quiero beber <u>jugo natural</u> de sandía. “I want to drink watermelon juice.”

Example c) **jugar** “*to play*” (links to **JUEGOS** “*games*,” **JUGUETES** “*toys*,” **DEPORTES** “*sports*”)

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Example d) **leer** “*to read*” (link to **LIBROS** “*books*”)

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3. Easy Access to Noun Categories in GRUPOS (word categories):

Mi bicicleta es azul. "My bicycle is blue."	hint: Find "bicicleta" in VEHÍCULOS or JUGUETES, and from that page COLORES for "es azul."
Quisiera ir a Hawaii. "I would like to go to Hawaii."	hint: Choose "quiero" and then you'll see "quisiera." Find "Hawaii" in GEOGRAFÍA/EE.UU.
Mi mejor amigo es mi perro. "My best friend is my dog."	hint: Choose "mejor" from DESCRIPCIONES, "amigo" from GENTE, and "perro" from ANIMALES.

4. Verb forms automatically agree with Pronouns selected:

Example a) Main page verb "comer" ("to eat")

Yo como. "I eat."
Tú comes. (<u>tú</u> is the informal form of "you") "You eat."
Usted come. (<u>usted</u> is the formal form of "you") "You eat." (<i>hint: find usted in GENTE</i>)
Él come. "He eats." (<i>hint: find él in GENTE</i>)
Nosotros comemos. "We eat." (<i>hint: find nosotros in GENTE</i>)
Ellos comen. "They eat." (<i>hint: find ellos in GENTE</i>)
Vosotros coméis. (<u>vosotros</u> is an informal form of plural "you," primarily used in Spain) "You eat." (<i>hint: find vosotros in GENTE</i>)
Note: **When you select ACCIONES and ACCIONES A-Z from the main page, the verb forms appear in the "present indicative" form (agrees with the <u>él/ella/usted</u> pronoun forms).

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Example c) Second page of ACCIONES “bailar” (“to dance”)

(hint: navigate using the arrow in bottom right of screen)

Yo bailo. “I dance.”
Tú bailas. “You dance.”
Él baila. “He dances.”
Nosotros bailamos. “We dance.”
Ellos bailan. “They dance.”

Example d) ACCIONES A-Z (hint: select the pronoun and then the “ACCIONES A-Z” link)

Yo conduzco. “I drive.”
Tú gastas un poco de dinero. (hint: find “dinero” in GRUPOS>DINERO) “You spend a lot of money.”
Nosotros peleamos. “We fight.”
Ellos disfrutan. “They enjoy.”

Notice there is an additional option with ACCIONES A-Z. In addition to choosing the pronoun on the main page and navigating to ACCIONES A-Z, once inside ACCIONES A-Z page you may select the pronoun on the left in yellow, and the verb forms change accordingly.

In your practice so far, you have probably noticed that the pronouns are not spoken when selected. In Spanish, subject pronouns typically are not necessary and are therefore often not spoken, as the verb form itself conveys the pronoun information. For example:

Hablo español.	(Yo is not necessary.)
I speak Spanish.	
Comemos carne.	(Nosotros is not necessary.)
We eat meat.	
Vives en México.	(Tú is not necessary.)
You live in Mexico.	

Using “WordPower48 Español**,” you can speak the pronoun for emphasis or clarity by selecting the underlined version of the pronoun that appears after the pronoun is selected.

5. Verb forms included in WordPower48 Español:

Present tense

Preterite	(past tense)
Imperfect	(was verbing)
Conditional	(would verb)
Future	(will verb)
Infinitive	(to + verb)
Present participle	(-ing)
Present/past perfect	(passive)

6. Verb forms appear in correct format automatically by context:

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Example c) Yo voy a ("I am going") + infinitives

(yo) Voy a (ACCIONES) leer. "I'm going to read."	(yo) Voy a jugar. "I'm going to play."	(yo) Voy a comer. "I'm going to eat."
---	---	--

(yo) Voy a (ACCIONES) comprarlo. "I'm going to buy it."	(comprar + lo + ←lo).
--	-----------------------

Example d) Yo he ("I have") + passive verb

(yo) He visto eso. "I have seen that."	*Note the words found after selecting buttons "este. . ." and "ese. . ." correspond to "this," "that," "these" and "those."
---	---

(yo) He regresado. "I have returned."	(yo) He terminado. "I have finished."	(yo) He ayudado "I have helped."
--	--	-------------------------------------

7. To access other verb tenses, select the pronoun and present tense of the verb. The preterite, imperfect, conditional and future tense forms appear on the top row.

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
8. Masculine/feminine nouns agree with Adjectives:

In Spanish, adjectives must match the noun in terms of number and gender. This is done automatically when choosing a color to modify a noun in categories such as **ANIMALES** (*animals*), **VEHÍCULOS** (*vehicles*), and **ROPA** (*clothes*). **Notice that the gender and number of the color changes automatically to agree with the nouns, and that the noun is positioned before the adjective.

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9. Plural fun!

In Spanish as in English, plural markers can be final –s or final –es depending upon the final letter of the noun. With two or three word noun concepts, the plural marker is typically on the first word. For example:

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11. Preguntas (Questions):

In Spanish, the proper punctuation for questions is an upside down question mark (¿) at the beginning of the sentence and a regular question mark (?) at the end. The blue (¿) capitalizes words such as que (what), dónde (where), por qué (why), cuánto (how much), etc., and the PREGUNTAS page remains **open**.

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12. Using the Imperative form of the Verb:

**Find imperative verbs by selecting “ACCIONES+IMPERATIVO” or “ACCIONES+INFINITIVO+IMPERATIVO”

Ven aquí. “Come here.”	(Select the “algo o todo pero” key on main page to find “aquí”)
Para ahora. “Stop now.”	(Find <u>ahora</u> in TIEMPO)
iComan! “You [plural] eat!”	(Select the “ellos/ellas/ustedes” key once you get to the IMPERATIVO page)

13. Other features in WordPower:

1. Select the GRUPOS key to find the CHISTES (Jokes). There are jokes and riddles. Try them out.
2. Select the GRUPOS key to find the NOTICIAS (News).
3. Select the GRUPOS key to find FOTOS. This gives access to a photo album that can be used to share pictures and experiences.

14. PUTTING IT ALL TOGETHER: try these sentences to practice your skills:

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