



We can pick our core word targets a variety of ways. With this form, we will be picking our core words and then the activities to go along with the targets. Literacy is always a great part of developing language so a book is the base here.



DECIDE to build communication
Determine the vocabulary
Explore Vocabulary on device
Change or alter as needed
Identify vocabulary with Button Capture in Chat Editor
Do practice
Expand with other activities

Start with a book _____

Determine the vocabulary. 80% Core/20% Fringe. Look for repeated words. Can the word that can be targeted in other activities <i>Brainstorm and list words here.</i>	Explore Vocabulary on Device. Change or alter as needed. Use another word in the system or add the word to the system <i>List the new word here.</i>	Identify the vocabulary path Use Button Capture in Chat Editor to create cheat sheet. WordPower 60.
Core: go <i>hunt</i> <i>catch</i> big one beautiful we not scared oh no!/uh-oh!/oops can't over under through it Fringe: bear, river, grass, mud	get I you cool	

<p>Words 80/20 rule List the target words identified in the reading activity</p>	<p>Activity/ Description Use the same words across other activities in the day.</p>	<p>Materials needed for additional activities (games, apps, etc)</p>	<p>What does the communication partner need to say and/or prompt to provide opportunity to communicate target vocabulary? Write a script or prompt reminder here.</p>	<p>Comments/ Data</p>
<p>Go Get Not Scared Under Over Through Cool Uh Oh! I You</p> <p>Bear River Grass Mud</p> 	<p>Activity 1: Practice the words -hunt for bear paw flash cards around the room</p>	<p>Sequences from Chat editor Bear paws Sticker</p>	<p>Ooh...we found one! Let's see where that word is!</p> <p>Model for them. Show them where to find it. Even act out the actions or describe words.</p>	
	<p>Activity 2: Human game-roll dice to see where they will go (grass/mud/river) and how they will get to the other side (under/over/through)</p>	<p>1 die- grass/mud/river 1 die- under/over/through Tall grass Towels for river Buckets for mud</p>	<p>Where will you go? Is that MUD or GRASS? Look like you will GO OVER! Where will you friend go? Let's put those words together. OVER GRASS That was cool!</p> <p>When someone else is going-you can tell him GO or Uh Oh!</p>	
	<p>Routine: Calendar time</p>	<p>Wall board Pointer</p>	<p>Let's say the days of the week- if it is NOT that day today, we can say NOT...Monday (NOT), Tuesday (NOT), Wednesday – GO!</p>	