






Unraveling the Mystery of AAC: Matching Systems to Needs Notes Handout

Trial Targets and Goals										
What are you measuring during the trial?	<table border="1" style="width: 100%; text-align: center;"> <tr> <td>Communicative functions</td> <td>Types of words used</td> <td># of words/symbols combined</td> </tr> <tr> <td>Prompting levels</td> <td>Access method</td> <td>Communication competencies</td> </tr> <tr> <td colspan="3">Other....?</td> </tr> </table>	Communicative functions	Types of words used	# of words/symbols combined	Prompting levels	Access method	Communication competencies	Other....?		
Communicative functions	Types of words used	# of words/symbols combined								
Prompting levels	Access method	Communication competencies								
Other....?										
What is the goal/priority for the AAC trial?	<table border="1" style="width: 100%; text-align: center;"> <tr> <td>Clarification or augmenting their verbal speech</td> <td>Expanding their utterances</td> <td>Increased initiation</td> </tr> <tr> <td>Moving beyond requesting</td> <td>Novel utterances</td> <td>Social exchanges</td> </tr> <tr> <td>Self-advocacy</td> <td>Increased independence</td> <td>Other...</td> </tr> </table>	Clarification or augmenting their verbal speech	Expanding their utterances	Increased initiation	Moving beyond requesting	Novel utterances	Social exchanges	Self-advocacy	Increased independence	Other...
Clarification or augmenting their verbal speech	Expanding their utterances	Increased initiation								
Moving beyond requesting	Novel utterances	Social exchanges								
Self-advocacy	Increased independence	Other...								
Where do you get this information?	<ul style="list-style-type: none"> Previous evaluations ALP for AAC Pragmatics Profile for AAC Users Discussions with team members <ul style="list-style-type: none"> SETT for AAC Participation Model Other sources you have available 									
What are you trialing on?	App version vs Trial Device									
Key Take Away	<ul style="list-style-type: none"> 									
Vocabulary File Selection										
Based on your assessment....	 <p>What access method were they successful with?</p> <p>How small of an item can they accurately select?</p> <p>What are they able to see?</p> <p>What is their language level?</p> <p>What words/vocabulary do they need access to right now?</p>									
Remember...	Access Before Communication! Never limit access to a robust vocabulary.									
Access Methods	<ul style="list-style-type: none"> Direct Select Keyguards or Touchguides Head Tracking (with or without headmouse) Eye Tracking Switch Access (1 or 2 switches) 									




Unraveling the Mystery of AAC:
Matching Systems to Needs Notes Handout

Keep in Mind	<ul style="list-style-type: none"> • Primary Message formulation • Consistent word placement • Amount of navigation to get to a word
Layout Size	<ul style="list-style-type: none"> • Accent Target Practice • PRC-Salttillo Size-wise Samples • Other:
Reflect larger vs smaller layout	
Strategies for Putting Access first	<ul style="list-style-type: none"> • Look at smaller button size first and step down the ladder if they need the bigger buttons • Think about the next 5 years, and the need for growing vocabulary • Motor automaticity • Try to create the least amount of change over time
Tools to Try	<ul style="list-style-type: none"> • Button borders • Hide/Mask buttons • Adding visual supports (e.g. magnifying, overlays, inverting colors) • Button timing • Keypadding/Button Margins • Explore touch tools • Look at the next step down in layouts
Key Take Aways	<ul style="list-style-type: none"> • •
Teach the Targets	
Starting point	Where is the student/client starting at and how are we getting them closer to the target?
Aided Language Stimulation	Speak AAC to Teach AAC!

Unraveling the Mystery of AAC:
Matching Systems to Needs Notes Handout

		<ul style="list-style-type: none"> • Speak it  • Say it with the device  • Wait  • Respond  	
AAC Prompting Hierarchy	<ul style="list-style-type: none"> • Start with the least amount of prompting and give an expectant pause. • Move to the next level if no response 		
	<ul style="list-style-type: none"> • Expectant Pause • Non-verbal Prompt • Environmental Cues • Ask a Question 	<ul style="list-style-type: none"> • Visual Prompt • Give Verbal Choices • Model on device/board • Partial Physical Prompt • Physical Prompt 	
Communication vs Compliance	<ul style="list-style-type: none"> • Communication, not compliance or testing • Communication should be fun, not drill and practice • Social communication before academic communication 		
	<p>“Child first has to learn to communicate before he can communicate to learn.” ~ Gail Van Tatenhove</p>		
	<p>Communication:</p> <ul style="list-style-type: none"> • What color do you like? • What animal do you want? • What would you like? More cookies? Or all done? 	<p>Compliance:</p> <ul style="list-style-type: none"> • What color is this? • What animal is this • Tell Me “more cookies, please” 	
Modeling without Expectation			
Modeling with Expectation with Wait Time			
Key Take Away	<ul style="list-style-type: none"> • 		
Data Collection and Comparison			

Unraveling the Mystery of AAC:
Matching Systems to Needs Notes Handout

Things to keep in mind...	<ul style="list-style-type: none"> • Does your funding source require a specific length of trial? • Does your funding source require you trial a specific number of devices? • Do you have enough data to show progress and make a recommendation
Baseline for data	<ul style="list-style-type: none"> • Know what level you are starting at for the targeted skill(s).
Data Collection	<ul style="list-style-type: none"> • Who?  • Where?  • How? 
Compare	<ul style="list-style-type: none"> • Look at your baseline data and post-trial data (after you taught the skill). • Look at progress made on each device/app trialed.
Key Take Aways	<ul style="list-style-type: none"> • •
Next Steps	<p>Keep trialing or ready to purchase?</p> <ul style="list-style-type: none"> • Did you try an iPad with an AAC app? • Did you do your due diligence for data collection? • Does your data support you making a recommendation?
	<p>Did you consider the least costly alternative?</p> <ul style="list-style-type: none"> • Low-tech • Mid-tech • Access to novel utterance generation • Independence • Functional communication needs
	<p>Make a decision as a team.</p>

References:
<ul style="list-style-type: none"> ▪ Aided Language Stimulation video by Chris Bugaj https://youtu.be/fIFNMky22-U ▪ ALP for AAC, https://www.alpforaac.com/ ▪ Byiers, B. J., Reichle, J., & Symons, F. J. (2012). Single-subject experimental design for evidence-based practice. <i>American Journal of Speech-Language Pathology</i>, 21(4), 397–414

Unraveling the Mystery of AAC:
Matching Systems to Needs Notes Handout

- Dada, S., & Alant, E. (2009). The effect of aided language stimulation on vocabulary acquisition in children with little or no functional speech. *American Journal of Speech-Language Pathology, 18*(1), 50–64.
- Finke, E. H., Davis, J. M., Benedict, M., Goga, L., Kelly, J., Palumbo, L., Peart, T., & Waters, S. (2017). Effects of a least-to-most prompting procedure on Multisymbol message production in children with autism spectrum disorder who use augmentative and alternative communication. *American Journal of Speech-Language Pathology, 26*(1), 81–98.
- Hyppa-Martin, J. (2016). AAC clinical education: Measuring outcomes with single subject designs. *Perspectives of the ASHA Special Interest Groups, 1*(12), 45–55.
- Korsten, J. (2002). Meaningful data - making sense of “+” and “-.” *Perspectives on Augmentative and Alternative Communication, 11*(3), 10–13.
- Thistle, J. J., Holmes, S. A., Horn, M. M., & Reum, A. M. (2018). Consistent symbol location affects motor learning in preschoolers without disabilities: Implications for designing augmentative and alternative communication displays. *American Journal of Speech-Language Pathology, 27*(3), 1010–1017.

Additional Notes: