


Chat Fusion 10



Introducing Chat Fusion, the newest communication device from Saltillo. Chat Fusion is driven by our powerful Chat software, but there is so much more. Fusion offers wireless charging capabilities, a front and rear camera, multiple access methods including optional headpointing and advanced scanning, and many other features.



KEY FEATURES

- Front and rear cameras
- SymbolStix or optional PCS symbols
- Word Prediction (Plus Models only)
- Extended battery life
- Wireless charging
- Single or dual switch scanning including advanced scanning features
- Chat Point (Optional headpointing)
- Ivona speech or optional Acapela speech

DEVICE SPECIFICATIONS

- Weight: 2.40 lbs.
- Size: 10.31" x 6.88" x 1.13"
- WUXGA capacitive touch screen
- 13 AHr Li-polymer rechargeable battery

VOICES AND ACCESSIBILITY

Accessibility:

- Finger touch
- Single or dual switch scanning

Voices:

- Headpointing (additional \$1295)
- Advanced scanning features
- IVONA synthesized speech
- Acapela speech (Additional \$195)

VOCABULARY

- | | | |
|----------------------------|------------------------|---------------------------|
| • WordPower 20 | • WordPower 48 Español | • MultiChat 15 Adolescent |
| • WordPower24 | • WordPower 60 | • MultiChat 15 Student |
| • WordPower24 with Phrases | • WordPower 80 | • MultiChat 15 Adult |
| • WordPower 30 | • WordPower 108 | • 4-Basic |
| • WordPower42 | • myCore | • Essence |
| • WordPower42 Basic | • myQuickChat | • Spelling |
| • WordPower 48 | • VocabPC™ | |



PRICING & ACCESSORIES

Device:

- Standard Model**
 #FUSION10-IVONA. \$4895
- Plus Model**
 #FUSION10-Plus-IVONA. \$5395
- Dedicated Plus Model** (includes 3 year warranty)
 #FUSION10-DPlus-IVONA. \$6395

Accessories:

- | | |
|---|--------------------------------------|
| Shoulder Strap \$25 | PCS Symbols \$375 |
| Carrying Case \$95 | Wheelchair mounting plate (included) |
| Keyguard (set of seven). \$175
(15, 20, 42, 48, 60, 80, 108 locations) | USB cable \$35 |
| Individual keyguards \$35 | Battery Charger \$69 |