






COMMUNICATION SOLUTIONS (see product flyers for more information)				
NOVA chat Line	Chat Fusion 10	TouchChat Express	TouchChat App	ChatWrap
5, 8, 10, 12 inches Dedicated Device	Dedicated Device Optional Chat Point	Mini & Standard Dedicated Device	Available through iTunes	Available on www.silver-kite.com
				

VOCABULARY OPTIONS

Vocabulary File Name	Type of Message Formulation				
	Spelling	Sentence	Phrase	-----Single Words-----	
				Nouns	Core
Communication Journey: Aphasia	Includes features & vocabulary to support people with aphasia – customization necessary				
Basic Scan		★	✓	✓	
Essence	★	✓	✓		
Spelling	★				
MyQuickChat	✓	★	✓	✓	
4 Basic			★	✓	
Vocab PC			★	✓	
MultiChat 15	✓	✓	★	★	✓
<input type="checkbox"/> Student <input type="checkbox"/> Adolescent <input type="checkbox"/> Adult					
WordPower: all configurations	✓	✓	✓	★	★
<input type="checkbox"/> 20 Simply <input type="checkbox"/> 24 <input type="checkbox"/> 24 with phrases <input type="checkbox"/> 30 <input type="checkbox"/> 42 Basic <input type="checkbox"/> 42 <input type="checkbox"/> 48 <input type="checkbox"/> 60 Basic <input type="checkbox"/> 60 <input type="checkbox"/> 80 <input type="checkbox"/> 108 <input type="checkbox"/> 108 with keyboard	-Configuration specific differences-		★		
			★		
				★	
			★		
myCore	✓	✓	✓	★	★

- ★: primary method of message formulation typically used
- ✓: additional methods of message formulation available

TouchChat includes additional vocabulary files available as in-app purchase
 Additional Language Options include: Spanish, Canadian French, Dutch, & German

ACCESS METHODS

Device	Direct Selection	Keyguards (standard & custom)	Scanning*	Headpointing with Chat Point
NOVA chat	X	X	X	
Chat Fusion	X	X	X	X
TouchChat Express	X	X	X (TouchChat scanning)	
TouchChat App	X	Available with ChatWrap	X (iOS scanning)	

Scanning software: switches (1, 2, screen), auditory prompts, visual effects

RESOURCES

CHAT EDITOR: FREE companion software for PC computer

How to get it

- Install the included CD or
- Download: <https://saltillo.com/products#chat-editor>

What to do with it

- Backup customized files
- Customize without taking the device away
- Create materials using Capture tool
- Learn vocabulary



ADDITIONAL RESOURCES

Make a **Low-Tech** Communication Book
<https://saltillo.com/chatcorner/content/29>

Attend a **Webinar**
<http://saltillo.com/webinars>

Watch Client and How To **Videos**
<http://youtube.com/saltillocorp>

Learn about **Realize Language & Data Logging**
www.realizelanguage.com

Listen to **Ivona & Acapela** Voices
www.ivona.com & www.acapela-group.com

Download Saltillo Pages for **News2You & Unique Learning**
<https://saltillo.com/curriculum-partners>
www.n2y.com

Learn about **TouchChat App**
<http://touchchatapp.com/>

Find **Implementation Ideas**
www.saltillo.com/chatcorner
www.practicalaac.org
www.inmaninnovations.com

Get **Support**:
 1-800-382-8622
<http://saltillo.com/support>

Join our **Newsletter**
<http://saltillo.com/newsletter>

Connect with Us   

NOTES