

Accessories and Pricing

	NOVA chat 12	NOVA chat 10	NOVA chat 8	NOVA chat 5
Standard Model	\$3995	\$3895	\$3495	\$2995
Plus Model	\$4495	\$4395	\$3995	\$3495
Dedicated Plus Model*	\$5495	\$5395	\$4995	\$4495

*Designed to meet the coverage requirements of the new NCD

Shoulder Strap \$25	Wheelchair mounting plate . . \$59 <small>(included with NOVA chat 12)</small>
Waist Belt \$25	USB Cable \$35
PCS Symbols \$375	Stylus \$10
Individual Keyguards . . . \$35	Battery Charger \$69
(See website for keyguard sets)	Carrying Case \$95



Family of Products

Support

Service:

Phone: 1-800-382-8622
Email: service@salttillo.com

Training:

<http://salttillo.com/training>

Implementation:

<http://salttillo.com/implementation>



2143 Township Rd 112
Millersburg, OH 44654
1.877.397.0178
info@salttillo.com
www.salttillo.com



FIND YOUR SALTILLO CONSULTANT

www.salttillo.com/rep



Living Life to its Fullest



www.salttillo.com

NOVA Chat Family


NOVA Chat is a family of Augmentative and Alternative Communication (AAC) devices designed by Salthillo Corporation. Built on an android platform, they offer the well known Chat software in a durable yet sleek portable design.



NOVA chat 12
Size: 12" x 8.38" x 1.25"
Weight: 2.95 lbs.



NOVA chat 10
Size: 9.75" x 7.13" x 1.63"
Weight: 1.95 lbs.



NOVA chat 8
Size: 8.5" x 5.13" x 1.63"
Weight: 1.4 lbs.




NOVA chat 5
Size: 5.5" x 2.69" x 1"
Weight: 0.55 lbs.

Chat Vocabulary

Chat devices ship with a number of default vocabulary files, or page sets. The vocabulary and organization of each vocabulary file was designed to address specific communication needs. All vocabulary files are fully customizable and can be personalized for each individual.


INCLUDED VOCABULARY FILES:

4 Basic	MultiChat 15 Student
MyCore	MultiChat 15 Adolescent
MyQuickChat	MultiChat 15 Adult
Essence	Vocab PC
Spelling	WordPower*



*WordPower

WordPower is a word-based vocabulary that is used to create sentences 'on the fly.' WordPower for NOVA chat was developed by Nancy L. Inman, M.A.T., CCC-SLP specifically for the NOVA chat.



WordPower is provided in a 20, 24, 30, 42, 48, 60, 80 and 108 button page layout.

CHAT EDITOR

Chat Editor is a supplemental program that has been provided for customizing the NOVA chat vocabulary on your Windows desktop or laptop computer.

Download now!
www.salthillo.com/products-chat-editor



KEY FEATURES

- Word Prediction (Plus Models only)
- Built-in front and rear cameras
- Single or dual switch scanning
- 4 rubberized color shields included
- SymbolStix or optional PCS symbols
- Ivona synthesized speech (including teen and child voices)