



TouchChat
Express

TouchChat™ Express

Introducing TouchChat™ Express, a new communication device from Saltillo. TouchChat™ Express combines the popular TouchChat app with the iPad and iPad mini. Now add industry leading service, support and training and you have the complete package.



KEY FEATURES

- TouchChat with WordPower
- Single or dual switch scanning
- Both English and Spanish included
- iShare (one year included)
- SymbolStix Symbols
- Word Prediction
- Acapela Speech
- Durable case with amplifier



• 1-877-397-0178 • www.saltillo.com •



DEVICE SPECIFICATIONS

iPad Air

- 6.88" x 9.63" x 1.68"
- 1.90 lbs

iPad mini

- 5.50" x 8.13" x 1.68"
- 1.35 lbs

VOICES AND ACCESSIBILITY

Accessibility:

- Finger touch
- Single or dual switch scanning

Voices:

- Acapela speech

VOCABULARY

English:

- WordPower
- VocabPC
- MultiChat 15
- MyCore
- Spelling
- 4 Basic
- MyQuickchat

Spanish:

- WordPower
- VocabPC
- MultiChat 15
- Teclado (Spelling)
- Sintaxis
- Sin Sintaxis



PRICING & ACCESSORIES

Device:

- EXPRESS \$3995
- EXPRESS-MINI . . . \$3995

Accessories:

Included with Purchase:

- Shoulder Strap
- Carrying Case
- Kick Stand

Available for Purchase:

- Individual keyguards \$35
(Keyguard sets available)
- Custom keyguard \$125
- Wheelchair mounting plate. . . \$59
- Waist Belt \$25