


Why would I want to print pages?

This offers a quick and easy way to create visual supports or low tech AAC options.

What pages might be good to print?

You can print any page of any vocabulary. Here are a few that you might want to consider



Home page of client's vocabulary. Name varies



Early Books Pages
WordPower 42 Basic, 60 Basic




TELL ME/MELD Boards in WordPower 60 Basic




Noun page with core (e.g., Art)
WordPower 42 Basic, 60 Basic, 108Keyboard

How do I print pages?



3. Check the pages you want to print.



Not sure what the page you want to print is called? Get into edit mode (if you are in an original, follow the prompts to quickly create a copy). Then you will see the name of any page at the top.

Download Chat Editor: <https://saltillo.com/products#chat-editor>