









Unraveling the Mystery of AAC:
Matching Systems to Needs Notes Handout

Feature Matching	
Options to Consider:	<ul style="list-style-type: none"> • Voice output or No Voice • Access Method <ul style="list-style-type: none"> ○ Direct Selection ○ Direct Selection with keyguard or touchguide ○ Switch scanning ○ Headtracking ○ Eyetracking • Primary Message Formulation: <ul style="list-style-type: none"> ○ Sentences ○ Phrases ○ Single Words/Core Words ○ Spelling ○ Specialty • Portability <ul style="list-style-type: none"> ○ Size of device ○ Weight of device • Durability • App vs Dedicated Device • Funding Sources <ul style="list-style-type: none"> ○ Out of pocket ○ Funding through insurance • Training and Support <ul style="list-style-type: none"> ○ During trial ○ After purchase
Key Take Aways	<ul style="list-style-type: none"> • •
PRC Solutions	
https://www.prentrom.com/products	
PRC Products	Accent Family of Devices (800, 1000, 1400) Prio Mini Via Pro Customer-bundled iOS device with LAMP Words for Life app
PRC Feature Match	<ul style="list-style-type: none"> • Access method- • Language System- • Durability- • Other Considerations-
	
PRC Trial device possibilities	   <p>Accent 1000 Accent 800 Via Pro with WFL Active</p>

Unraveling the Mystery of AAC:
Matching Systems to Needs Notes Handout

<p>Unity Points</p>	<ul style="list-style-type: none"> • Unity configurations: <ul style="list-style-type: none"> • 4, 8, 15, 28, 36, 45, 60, 84 and 144 locations • 1-hit • Sequenced • Vocabulary is learned through motor planning and consistent patterns • Word finder • Vocabulary builder
<p>PRC notes</p>	
<p>Salttillo Solutions</p>	
<p>https://salttillo.com/products</p>	
<p>Salttillo Products</p>	<p>NovaChat Family of devices (5, 8, 10, 12)</p> <p>TouchChat Express 8</p> <p>Via Pro with TouchChat</p> <p>Customer-bundled iOS device with TouchChat app</p>
<p>Salttillo Feature Match</p> 	<ul style="list-style-type: none"> • Access Method- • Language System- • Durability- • Other Considerations-
<p>Salttillo Trial Device possibilities</p>	<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;">  <p>NovaChat 8 Active</p> </div> <div style="text-align: center;">  <p>NovaChat 10 Active</p> </div> <div style="text-align: center;">  <p>Via Pro with TouchChat Active</p> </div> </div>
<p>WordPower Key Points</p>	<ul style="list-style-type: none"> • WordPower configurations: <ul style="list-style-type: none"> ○ 20, 25, 42, 48, 60, 80, 108 and 140 locations • Vocabulary uses high frequency core words for creating novel utterances. • Logical next word navigation • Word finder
<p>Salttillo Notes</p>	

Unraveling the Mystery of AAC:
Matching Systems to Needs Notes Handout

Next Steps	
Contact regional rep's to find out about:	<ul style="list-style-type: none"> Trial devices Funding process Training and support during the trial Training and support after purchase
You can start at Explore AAC for your PRC and Saltillo solutions- https://exploreaac.com/	
PRC Consultant	<ul style="list-style-type: none">
Saltillo Consultant	<ul style="list-style-type: none">

References:
<ul style="list-style-type: none"> Explore AAC, https://exploreaac.com/ PRC Product Solutions, https://www.prentrom.com/products Saltillo Product Solutions, https://saltillo.com/products
Additional Notes: