

Name: _____

Date: _____

Access Method	Direct Selection	Keyguard Visual supports	Touchguide Touch settings (dwell, delay, release)	stylus
	Mouse	Head mouse Joystick Trackball	light pointer trackpad alternative	
	Scanning	Method *1-switch, 2 -switch, auto) Feedback (outline, magnify, invert colors, auditory, combo) Pattern (linear, row/column, block, flow, combo) Switch type (button, proximity, voice, grip, other)		
	Head Tracking	With head mouse With facial recognition		
	Eye Tracking	Interaction settings (calibration/camera settings) Activation signal (dwell, blink, switch) Feedback (color, cursor, style, outline, highlight)		
Primary Message Formulation	Sentence based Single/Core Words			
Size	_____ inch screen			
Durability	Yes _____ No _____			
Portability	handle	carrying strap	carrying case	
Positioning	kickstand	table mount	wheelchair mount	
Integration	Environmental controls Internal computer access External computer access MP3 player Calculator			
Dedicated	Dedicated device		App-based	
Funding Source	Private insurance _____ MA/Medicare _____ Out-of-pocket			
Features to Prioritize:				
Notes:				

Trial Options

Device 1	Device Name:	
	Vocabulary Set:	
	Accessories Needed:	
	Rep contact info:	
Device 2	Device Name:	
	Vocabulary File:	
	Accessories Needed:	
	Rep contact info:	
Device 3	Device Name:	
	Vocabulary File:	
	Accessories Needed:	
	Rep contact info:	

Name: Max C.

Date: 1-28-21

Access Method	Direct Selection	<u>Keyguard</u>	<u>Touchguide</u>	stylus
		Visual supports	<u>Touch settings (dwell, delay, release)</u>	
	Mouse	Head mouse Joystick Trackball	<u>light pointer</u> <u>trackpad</u> <u>alternative</u>	
	Scanning	Method *1-switch, 2-switch, auto) Feedback (outline, magnify, invert colors, auditory, combo) Pattern (linear, row/column, block, flow, combo) Switch type (button, proximity, voice, grip, other)		
	Head Tracking	With head mouse With facial recognition		
	Eye Tracking	Interaction settings (calibration/camera settings) Activation signal (dwell, blink, switch) Feedback (color, cursor, style, outline, highlight)		
Primary Message Formulation	<u>Sentence based</u> <u>Single/Core Words</u>	<u>Phrase Based</u> <u>Spelling</u>		
Size	<u>8-10</u> inch screen			
Durability	<u>Yes</u>	No		
Portability	<u>handle</u>	<u>carrying strap</u>	carrying case	
Positioning	<u>kickstand</u>	table mount	wheelchair mount	
Integration	Environmental controls		Internal computer access	
	External computer access			
	<u>MP3 player</u>	Calculator		
Dedicated	<u>Dedicated device</u>	App-based		
Funding Source	Private insurance <u>BCBS</u>			
	MA/Medicare <u>WI Medicaid</u>			
	Out-of-pocket			
Features to Prioritize:	<ul style="list-style-type: none"> • durable = low tone, falls/trips, bumps into things • Core words = has used manual core boards well, need novel utterances • direct select with support = needs help for accuracy 			
Notes:	<ul style="list-style-type: none"> - uses an iPad at school, but want to see if smaller size would be more portable too - higher receptive lang. so we want enough vocab to grow with him 			



Device	Accent® Family of Devices	PRiO® Mini	Via® Pro	Customer-Bundled iOS Device
Operating System	Windows	Apple iOS	Apple iOS	Apple iOS
Hardware	<ul style="list-style-type: none"> - Splash & dust resistant design - Exclusively manufactured tablet - Chemically hardened glass - Impact case w/magnesium frame - Available in 8," 10," 14" sizes - Xtreme™ case available for added protection (Accent 1000 & 800) 	<ul style="list-style-type: none"> - 16 GB iPad Mini - Protective classic wrap case with built-in Bluetooth speakers - Dual charger and multiple switch ports 	<ul style="list-style-type: none"> - 256 GB, 11" iPad Pro - Protective Active or Classic case with built-in Bluetooth speakers - Single charger and multiple switch ports - Splash & dust resistant design - Sturdy handle and adjustable stand 	<ul style="list-style-type: none"> - Any iPad® or iPad Mini™ - Consumer available cases and speakers
Vocabulary Options	<ul style="list-style-type: none"> - Unity®, LAMP Words for Life® app, and CoreScanner™ - UNIDAD® - Essence® & WordPower™ app - Engage™ & Page-Based 	<ul style="list-style-type: none"> - LAMP Words for Life® app 	<ul style="list-style-type: none"> - LAMP Words For Life® app - TouchChat® HD with WordPower™ app - In-app vocabulary options available (not available for all funding sources) 	<ul style="list-style-type: none"> - LAMP Words for Life® app - WordPower™ app
Product Training/Support	<ul style="list-style-type: none"> - Free support from AAC Consultants - Free 1-year subscription to Realize Language & AAC Language Lab - Live webinars - On-demand trainings - YouTube tutorials 	<ul style="list-style-type: none"> - Free support from AAC Consultants - Free 1-year subscription to Realize Language & AAC Language Lab - Live webinars - On-demand trainings - YouTube tutorials 	<ul style="list-style-type: none"> - Free support from AAC Consultants - Free 1-year subscription to Realize Language & AAC Language Lab - Live webinars - On-demand trainings - YouTube tutorials 	<ul style="list-style-type: none"> - Fee-based support from AAC Consultant - Fee-based subscription to Realize Language and AAC Language Lab - Live webinars - On-demand trainings - YouTube tutorials
Technical Support	<ul style="list-style-type: none"> - Free phone and online support 	<ul style="list-style-type: none"> - Free phone and online support 	<ul style="list-style-type: none"> - Free phone and online support 	<ul style="list-style-type: none"> - Email support - Fee-based premium support plan
Alternative Access	<ul style="list-style-type: none"> - Keyguards and TouchGuides - Switch scanning (wired & wireless) - Head Tracking - Eye Tracking 	<ul style="list-style-type: none"> - 84 location keyguard - Switch scanning (wired & wireless) 	<ul style="list-style-type: none"> - 84 location keyguard & TouchGuide - Switch scanning (wired & wireless) - Head Tracking (within TouchChat) 	<ul style="list-style-type: none"> - Consumer available keyguards - Bluetooth switch only
Warranty	<ul style="list-style-type: none"> - 2 years 	<ul style="list-style-type: none"> - 1 year 	<ul style="list-style-type: none"> - 1-2 years 	<ul style="list-style-type: none"> - N/A
Price	<ul style="list-style-type: none"> - \$6,595 to \$7,595, varies by device 	<ul style="list-style-type: none"> - \$3,995 	<ul style="list-style-type: none"> - \$4,595 - \$5,495, varies by device 	<ul style="list-style-type: none"> - Varies, but typically between \$600-\$1,200
Fundability	<ul style="list-style-type: none"> - Insurance, Medicaid, Medicare Funding 	<ul style="list-style-type: none"> - Insurance, Medicaid, Medicare Funding 	<ul style="list-style-type: none"> - Insurance, Medicaid, Medicare Funding 	<ul style="list-style-type: none"> - Not fundable

Additional Information Find your PRC Regional Consultant at prentrom.com/consultants, Training for Accent Products/Programming, Unity & LAMP Words for Life Language Systems – prentrom.com/training, Information on the LAMP Therapy Approach from The Center for AAC and Autism – aacandautism.com

800-262-1984
prentrom.com



























































Device	NovaChat® Family of Devices	TouchChat® Express 8	Via® Pro	Customer-Bundled iOS Device
Operating System	Android	Apple iOS	Apple iOS	Apple iOS
Hardware	<ul style="list-style-type: none"> - Splash & dust resistant design - Protective Active or Classic wrap case with built-in Bluetooth speakers - Sturdy handle & adjustable stand - Available in 5," 8," 10," 12" sizes 	<ul style="list-style-type: none"> - 16 GB iPad Mini - Protective classic wrap case with built-in Bluetooth speakers - Dual charger and multiple switch ports 	<ul style="list-style-type: none"> - 256 GB, 11" iPad Pro - Protective Active or Classic case with built-in Bluetooth speakers - Single charger and multiple switch ports - Splash & dust resistant design - Sturdy handle and adjustable stand 	<ul style="list-style-type: none"> - iPad®, iPad mini™, iPhone®, and iPod touch® - Consumer available cases and speakers
Vocabulary Options	<ul style="list-style-type: none"> - All standard Saltillo vocabulary options, including the full WordPower™ Suite 	<ul style="list-style-type: none"> - TouchChat HD with WordPower™ app - In-app vocabulary options available (not available for all funding sources) 	<ul style="list-style-type: none"> - TouchChat® HD with WordPower™ app - In-app vocabulary options available (not available for all funding sources) - LAMP Words For Life® app 	<ul style="list-style-type: none"> - TouchChat® HD with WordPower™ app - In-app vocabulary options available (not available for all funding sources)
Product Training/ Support	<ul style="list-style-type: none"> - Free support from AAC Consultants - Free 1-year subscription to Realize Language & AAC Language Lab - Live webinars - On-demand trainings - YouTube tutorials 	<ul style="list-style-type: none"> - Free support from AAC Consultants - Free 1-year subscription to Realize Language & AAC Language Lab - Live webinars - On-demand trainings - YouTube tutorials 	<ul style="list-style-type: none"> - Free support from AAC Consultants - Free 1-year subscription to Realize Language & AAC Language Lab - Live webinars - On-demand trainings - YouTube tutorials 	<ul style="list-style-type: none"> - Fee-based support from AAC Consultant - Fee-based subscription to Realize Language and AAC Language Lab - Live webinars - On-demand trainings - YouTube tutorials
Technical Support	<ul style="list-style-type: none"> - Free phone and online support 	<ul style="list-style-type: none"> - Free phone and online support 	<ul style="list-style-type: none"> - Free phone and online support 	<ul style="list-style-type: none"> - Email support - Fee-based premium support plan
Alternative Access	<ul style="list-style-type: none"> - Keyguards and TouchGuides - Switch Scanning (wired & wireless) 	<ul style="list-style-type: none"> - Keyguards and TouchGuides - Switch Scanning (wired & wireless) 	<ul style="list-style-type: none"> - Keyguards and TouchGuides - Switch Scanning (wired & wireless) - Head Tracking (within TouchChat) 	<ul style="list-style-type: none"> - Consumer available keyguards - Bluetooth switch only
Warranty	<ul style="list-style-type: none"> - 1-3 years 	<ul style="list-style-type: none"> - 1 year 	<ul style="list-style-type: none"> - 1-2 years 	<ul style="list-style-type: none"> - N/A
Price	<ul style="list-style-type: none"> - \$4,595 - \$5,495, varies by device 	<ul style="list-style-type: none"> - \$3,995 	<ul style="list-style-type: none"> - \$4,595 - \$5,495, varies by device 	<ul style="list-style-type: none"> - Varies, but typically between \$600-\$1,200
Fundability	<ul style="list-style-type: none"> - Insurance, Medicaid, Medicare Funding 	<ul style="list-style-type: none"> - Insurance, Medicaid, Medicare Funding 	<ul style="list-style-type: none"> - Insurance, Medicaid, Medicare Funding 	<ul style="list-style-type: none"> - Not fundable



what 	when 	where 	ABC 123	PEOPLE 	QUESTN 	ACTIONS 	SOCIAL 	PLACES 	TIME 	GROUPS 	DESCRB 
I	me 	how 	who 	why 	again 	please 	thank you 	problem 	now 	bad 	good 
my/mine 	am	to	be	feel 	give 	listen 	happy 	sad 	tired 	okay 	cool 
it 	is are	will	come 	hurt 	hear 	know 	that 	a	the	and 	more 
you 	can	eat 	drink 	finish 	get 	love 	make 	need 	all 	at	some 
your 	do	go 	help 	open 	put 	say/talk 	see/look 	first 	then 	for of	on 
here 	have 	like 	play 	read 	stop 	walk 	show 	wait min 	in 	up 	off 
yes 	no/don't 	want 	take 	tell 	turn 	watch 	wear 	work 	out 	down 	with

	Hi 	What's up? 	You? 	Fine. 	Not good. 	Nice meet 	I'm going 	Goodbye 	
GIVE 	NAMES 	SOCIAL 	DETERMINER 	CONJUNCTIO. 	PREPOSITI... 	QUESTIONS ? word	a	+s	AGAIN 
I 	WE 	are 	is 	were 	was 	to 	KEYBOARD 	an 	the 
YOU 	THEY 	PLAY 	LIKE 	WORK 	HAVE 	FEEL 	READ 	MORE 	STOP 
IT 	HE 	WANT 	COME 	GOOD 	DO 	LIVE 	GO 	COLOR 	HELP 
CLEAR 	SHE 	CHANGE 	EAT 	MAKE 	BAD 	DRINK 	TURN 	PUT 	SLEEP 